

# Maximum oxygen transfer capacity



AppNote by **Kuhner** shaker

Keebler, M.V.<sup>1</sup>, Laidlaw, R.D.<sup>1</sup>, Anderlei, T.A.<sup>2</sup>

1. Kuhner Shaker, Inc., San Carlos, CA, United States
2. AG Kuhner, Birsfelden, Switzerland

Oxygen is an essential nutrient for all aerobic organisms and knowing the maximum oxygen transfer capacity ( $OTR_{max}$ ) of your cultivation system is critical for experimental planning and optimization. Here we discuss the key parameters that define the  $OTR_{max}$  and illustrate how these parameters can be used to theoretically estimate the  $OTR_{max}$  for any aerobic system. We also discuss the use of proportionality factors to relate empirical measurements of the  $OTR_{max}$  from a sulfite system to any biological system.

## What is the $OTR_{max}$ and why is it important?

The  $OTR_{max}$  defines the maximum capacity for oxygen transfer between the gas and the liquid inside a shaking vessel [1]. Once gaseous oxygen has dissolved into the culture media, it is available as a nutrient for living organisms to consume. The rate at which dissolved oxygen is consumed in the liquid by aerobic organisms is referred to as the oxygen uptake rate (OUR). When the OUR is greater than the  $OTR_{max}$ , the metabolic demand for oxygen is not being met and growth of the culture will be limited by the  $OTR_{max}$  [1-5]. This phenomenon, known as 'oxygen limitation', can be addressed by optimizing the shaking conditions to ensure that the  $OTR_{max}$  is greater than the OUR requirements of the culture.

## How is the $OTR_{max}$ calculated?

The key parameters that will impact the  $OTR_{max}$  of any shaken process include the: (1) osmolality of the solution, (2) shaking frequency, (3) filling volume, (4) shaker diameter, and (5) vessel diameter. Meier et al. (2016) describe a model for the calculation of  $OTR_{max}$  in shake flasks, utilizing the equation below with the above mentioned parameters [5]. The equation also includes a variable for the absolute pressure in the flask ( $p_R$ ), which is assumed to be equivalent to ambient pressure (1 bar), and the mole fraction of  $O_2$  in the ambient air ( $y_{O_2}$ ), which is assumed to be ~0.21. This equation can be utilized as a valuable tool to calculate and compare the  $OTR_{max}$  values for different media under the same shaking conditions. For instance, if a user were growing Chinese Hamster Ovarian (CHO) cells in 40mL of ProCHO™ 5 media, by shaking in a 250mL Erlenmeyer (non-baffled) flask at 150rpm with a 50mm orbital diameter, they should expect an  $OTR_{max}$  value of ~10.8 mmol/L/h.

If the user then decided to switch to the ExpiCHO™ Stable Production Medium and maintain the same shaking conditions, we would only need to change the osmolality value in the equation from 0.350 Osmol/kg [6] to 0.295 Osmol/kg [7]. This would increase the  $OTR_{max}$  to ~11.0 mmol/L/hr. By adjusting the values for the different parameters in this equation, the user can easily predict the  $OTR_{max}$  across different process conditions. Table 1 lists a few of these calculations for different media. These predictions can also be empirically confirmed using an online, off-gas analyzer like the Kuhner TOM ([www.kuhner.com](http://www.kuhner.com)).

## Converting $OTR_{max}$ between systems

The sodium sulfite system is a useful tool for mimicking biological activity in a shaking vessel [4]. The chemical reaction between sulfite ( $SO_3$ ) and dissolved oxygen ( $O_2$ ) is catalyzed with cobalt to produce sulfate ( $SO_4$ ) and consume oxygen. These chemical systems have been widely used to study how the  $OTR_{max}$  varies with respect to different shaking conditions [4,5,8,9]. Many users still rely on these data sets for bioprocess development and optimization, and the  $OTR_{max}$  values from these chemical systems can be easily converted to any biological media using a proportionality factor. Since the equation listed below is specific to Erlenmeyer flasks, this approach of relating empirical measurements from a chemical system to a biological system could be useful for other vessel types or to capture additional variables that may not be factored into the equation.

These proportionality factors usually range between 1 to 3 and will likely vary with media composition and shaking conditions (Figure 1). Table 2 (Appendix) lists a few examples of proportionality factors for different media and process conditions.

$$OTR_{max} = 3.72 \times 10^{-4} \times \text{Osmol}^{0.05} \times n^{(1.18 - \frac{\text{Osmol}}{10.1})} \times V_L^{-0.74} \times d_0^{0.33} \times d^{1.88} \times p_R \times y_{O_2}$$

|   |                          |                      |                      |                    |
|---|--------------------------|----------------------|----------------------|--------------------|
| Osmolality<br>(Osmol kg <sup>-1</sup> ) | Shake Frequency<br>(RPM) | Filling Vol.<br>(mL) | Shaking Dia.<br>(cm) | Flask Dia.<br>(mm) |
|---|--------------------------|----------------------|----------------------|--------------------|

# Maximum oxygen transfer capacity

AppNote by **Kuhner** shaker




**Table 1. Calculated  $OTR_{max}$  values for different media.**

$OTR_{max}$  calculations for different media osmolality and shaking conditions, based on equation described in Meier et al., (2016) [5].

| Media                             | $OTR_{max}$ (mmol/L/h) | Constant | Osmolality (Osmol/kg) | Shake Frequency (RPM) | Filling Volume (mL) | Shaker Diameter (cm) <sup>a</sup> | Flask Diameter (mm) | Absolute Pressure (bar) | O <sub>2</sub> mole fraction |
|-----------------------------------|------------------------|----------|-----------------------|-----------------------|---------------------|-----------------------------------|---------------------|-------------------------|------------------------------|
| ProCHO 5 Media                    | 10.8                   | 3.72E-04 | 0.350                 | 150                   | 40                  | 5                                 | 85                  | 1                       | 0.21                         |
| ExpiCHO™ Stable Production Medium | 11.1                   | 3.72E-04 | 0.295                 | 150                   | 40                  | 5                                 | 85                  | 1                       | 0.21                         |
| ExpiCHO™ Expression Medium        | 11.0                   | 3.72E-04 | 0.295                 | 150                   | 40                  | 5                                 | 85                  | 1                       | 0.21                         |
| Expi293™ Expression Medium        | 11.0                   | 3.72E-04 | 0.285                 | 150                   | 40                  | 5                                 | 85                  | 1                       | 0.21                         |
| TB-Medium                         | 62.5                   | 3.72E-04 | 0.500                 | 300                   | 10                  | 5                                 | 85                  | 1                       | 0.21                         |
| LB-Medium                         | 69.8                   | 3.72E-04 | 0.240                 | 300                   | 10                  | 5                                 | 85                  | 1                       | 0.21                         |
| 1M Sodium Sulfite System          | 24.4                   | 3.72E-04 | 2.300                 | 300                   | 10                  | 5                                 | 85                  | 1                       | 0.21                         |
| 0.5M Sodium Sulfite System        | 41.7                   | 3.72E-04 | 1.300                 | 300                   | 10                  | 5                                 | 85                  | 1                       | 0.21                         |

<sup>a</sup> Average width of 250mL, unbaffled flask.



**Figure 1. Proportionality of  $OTR_{max}$  between culture media and sodium sulfite systems.**

Data represent ratio of the  $OTR_{max}$  for *E. coli* in terrific broth (TB) media (squares) or Wilms-MOPS (WM) media (circles) with the  $OTR_{max}$  measured with a 1M sulfite system (blue) or a 0.5M sulfite system (red) as a function of shaking frequency. Data used with permission from Meier et al., (2016) [5].

In summary, the  $OTR_{max}$  defines the capacity for oxygen transfer in a system and is an important parameter for the scaling or optimization of any aerobic process. The model described by Meier et al., 2016 [5] can be used to calculate the theoretical  $OTR_{max}$  of any system and this can be empirically validated using an online, off-gas analyzer like the Kuhner TOM ([www.kuhner.com](http://www.kuhner.com)). Many  $OTR_{max}$  measurements have come from non-biological systems and are often used to compare the performance of different shaking conditions. These  $OTR_{max}$  measurements can be easily translated into any biological system using a proportionality factor for that system. Using these tools to determine the  $OTR_{max}$  of a system will help ensure that the OUR of a system is never greater than the  $OTR_{max}$ .

## References:

- Hirose, Y., Sonoda, H., Kinoshita, K., Okada, H. 1965. Studies on oxygen transfer in submerged fermentations part IV. Determine of oxygen transfer rate and respiration rate in shaken cultures using oxygen analyzers. *Agr. Biol. Chem.* 30(1): 49-58.
- Van Suijdam, J.C., Kossen, N.W.F., Joha, A.C. 1978. Model for oxygen transfer in a shake flask. 1978. *Biotechnol. Bioeng.* 20(11): 1695-1709.
- Büchs, J. 2001. Introduction to advantages and problems of shaken cultures. *Biochem. Eng. J.* 7 (2): 91-98.
- Maier, U., Losen, M., Büchs, J. 2004. Advances in understanding and modeling the gas-liquid mass transfer in shake flasks. *Biochem. Eng. J.* 17(3): 155-167.
- Meier, K., Klockner, W., Bonhage, B., Antonov, E., Regestein, L., Büchs, J. 2016. Correlation for the maximum oxygen transfer capacity in shake flasks for a wide range of operating conditions and for different culture media. *J. Biosci. Bioeng.* 109: 228-235.
- Lonza Verviers, S.p.r.l. 2007. BioWhittaker™ ProCHO™ 4 and ProCHO™ 5 powder kits. Document: BE005-1.
- Gibco. 2020. Cat. #: A37110. Certificate of Analysis. Lot #: 2216244.
- Kensy, F., Zimmermann, H.F., Knabben, I., Anderlei, T., Trauthwein, H., Dingerdissen, U., Büchs, J. 2005. Oxygen transfer phenomena in 48-well microtiter plates: determination by optical monitoring of sulfite oxidation and verification by real-time measurement during microbial growth. *Biotechnol. Bioeng.* 89(6): 698-708.
- Hermann, R., Walther, N., Maier, U., Büchs, J. 2001. Optical method for the determination of the oxygen-transfer capacity of small bioreactors based on sulfite oxidation. *Biotechnol. Bioeng.* 74(5): 355-63.
- Giese, H., Kruthof, P., Meier, K., Sieben, M., Antonov, E., Hommes, R.W., Büchs J. 2014. Improvement and scale-down of a *Trichoderma reesei* shake flask protocol to microtiter plates enables high-throughput screening. *Biotechnol. Bioeng.* 118(6): 702-9.
- Seletzky, J.M., Noak, U., Fricke, J., Welk, E., Eberhard, W., Knocke, C., Büchs, J. 2007. Scale-up from shake flasks to fermenters in batch and continuous mode with *Corynebacterium glutamicum* on lactic acid based on oxygen transfer and pH. *Biotech. Bioeng.* 98(4): 800-11.
- Stöckmann C, Maier U, Anderlei T, Knocke C, Gellissen G, Büchs J. 2003. The oxygen transfer rate as key parameter for the characterization of *Hansenula polymorpha* screening cultures. *J. Ind. Microbiol. Biotechnol.* 30(10): 613-622.
- Rahmen, N., Fulton, A., Ihling, N., Magni, M., Jaeger, K.E., Büchs, J. 2015. Exchange of single amino acids at different positions of a recombinant protein affects metabolic burden in *Escherichia coli*. *Microb. Cell Fact.* 14:10.
- Scheidle, M., Dittrich, B., Klinger, J., Ikeda, H., Klee, D., Büchs, J. 2011. Controlling pH in shake flasks using polymer-based controlled-release discs with pre-determined release kinetics. *BMC Biotechnol.* 11:25.
- Stöckmann C., Palmen, T.G., Schroer, K., Kunze, G. Gellissen, G., Büchs, J. 2014. Definition of culture conditions for *Arxula adenivorans*, a rational basis for studying heterologous gene expression in this dimorphic yeast. *J. Ind. Microbiol. Biotechnol.* 41(6): 965-76.
- Wilming, A., Begemann, J., Kuhne, S., Regestein, L., Bongaerts, J., Evers, S., Maurer, K.H., Büchs, J. 2013. Metabolic studies of  $\gamma$ -polyglutamic acid production in *Bacillus licheniformis* by small-scale continuous cultivations. *J. Ferm. Bioeng.* 73: 29-37.
- Stöckmann, C., Losen, M., Dahlems, U., Knocke, C., Gellissen, G., Büchs, J. 2003. Effect of oxygen supply on passaging, stabilizing and screening of recombinant *Hansenula polymorpha* production strains in test tube cultures. *FEMS Yeast Res.* 4(2): 195-205.
- Duetz, W.A., Rüedi, L., Hermann, R., O'Connor, K., Büchs, J., Witholt, B. 2000. Methods for intense aeration, growth, storage, and replication of bacterial strains in microtiter plates. *Appl. Environ. Microbiol.* 66(6): 2641-6.
- Giese, H., Azizan, A., Kümmel, A., Liao, A., Peter, C.P., Fonseca, J.A., Hermann, R., Duarte, T.M., Büchs, J. 2014. Liquid films on shake flask walls explain increasing maximum oxygen transfer capacities with elevating viscosity. *Biotechnol. Bioeng.* 111(2): 295-308.

# Maximum oxygen transfer capacity

AppNote by **Kuhner** shaker



## Appendix

**Table 2. Proportionality factors to convert the OTR<sub>max</sub> measurements from a sulfite system to a biological media.**

Proportionality factors relating OTR<sub>max</sub> measurements between different culture media and different sodium sulfite systems are listed along with the associated culture parameters.

| [Sulfite] | Orbit (mm) | Shaking Frequency (RPM) | Shaking Vessel                   | Filling Volume (%) | Cell Type                     | Media                | Proportionality Factor | Literature |
|-----------|------------|-------------------------|----------------------------------|--------------------|-------------------------------|----------------------|------------------------|------------|
| 1M        | 50mm       | 50-500rpm               | Shake flask (250mL)              | 4-16               | <i>P. pastoris</i>            | YPG                  | 2.8                    | [4]        |
| 1M        | 50mm       | 200-350                 | Shake flask (250mL), MTP (24w)   | 4-20, 21-29        | <i>T. reesei</i>              | Complex media        | 1.9 to 2.4             | [10]       |
| 1M        | 50mm       | 300                     | Shake flask (500mL)              | 8                  | <i>C. glutamicum</i>          | Complex media        | 1.8                    | [11]       |
| 1M        | 50mm       | 300                     | Shake flask (500mL)              | 8                  | <i>H. polymorpha</i>          | YNB                  | ~1.52                  | [12, 11]   |
| 1M        | 50mm       | 350                     | Shake flask (250mL)              | 4                  | <i>E. coli</i> (BL21-DE3)     | Terrific Broth       | ~1.97 <sup>a</sup>     | [13, 10]   |
| 1M        | 50mm       | 350                     | Shake flask (250mL)              | 4                  | <i>E. coli</i> (BL21-DE3)     | Wilms-MOPS           | ~1.59 <sup>a</sup>     | [13, 10]   |
| 1M        | 50mm       | 350                     | Shake flask (250mL)              | 4                  | <i>E. coli</i> (K12 and BL21) | Wilms-MOPS           | ~1.47 <sup>a</sup>     | [14, 10]   |
| 1M        | 50mm       | 350                     | Shake flask (250mL)              | 4                  | <i>A. adenivorans</i> (LS3)   | SYN6-PO <sub>4</sub> | ~1.43 <sup>a</sup>     | [15, 10]   |
| 0.5M      | 25mm; 50mm | 450; 250, 300           | MTP (48w)                        | 29                 | <i>H. polymorpha</i>          | YNB                  | ~1.35 <sup>a</sup>     | [8]        |
| 0.5M      | 50mm       | 350                     | Shake flask (250mL)              | 4                  | <i>B. licheniformis</i>       | V3 MOPS              | 1.84                   | [16]       |
| 0.5M      | 50mm       | 320 and 240             | Tube (15mL), Shake flask (250mL) | 6, 4               | <i>H. polymorpha</i>          | YNB                  | 1.5                    | [17]       |
| 0.5M      | 50mm       | 300                     | MTP (96w) w/System Duetz         | 21                 | <i>P. putida</i>              | Mineral medium       | ~1.38                  | [18]       |
| 0.5M      | 25mm       | 250                     | Shake flask (250mL)              | 12                 | <i>P. putida</i>              | Mineral medium       | ~1                     | [9]        |
| 0.35M     | 50mm       | 350                     | Shake flask (250mL)              | 12                 | <i>C. glutamicum</i> (DM1730) | Complex media        | ~1.1                   | [19]       |

<sup>a</sup> OTR<sub>max</sub> values utilized for these proportionality factors were separately reported, though measured under the same conditions