




CO₂ stability inside the vessel

What happens when the Beluga door opens?



ISF1-Z Basic Beluga
Perfect fit for automation



Why this matters

Researchers often worry that frequent door openings in CO₂-controlled incubator shakers will disturb pH-relevant CO₂ levels inside their cultures. While chamber CO₂ can change rapidly, the key question is: **What happens inside the cultivation vessel itself?**

Measurement approach

Using the **Kühner TOM** off-gas analysis system, CO₂ was measured directly in the headspace of **600 mL TPP TubeSpin® bioreactors** during shaking via a bypass—without vessel-internal sensors or wires.

What actually happens when the door opens

Opening the Beluga door for 30 seconds cut the chamber CO₂ from 8% to ~4%. Inside the tubes, however, CO₂ decreased only ~10%, stabilizing around 7.3% CO₂. The response was **slow, modest, and nearly identical** across all vessel positions (B). Shorter openings (15 s) caused even smaller changes (C).

Repeated door openings

Frequent openings led to a **stable oscillating equilibrium** inside the tubes (D+E+F)

- 3 - minute breaks → ~6.4% CO₂
- 8 - minute breaks → ~7.3% CO₂
- 15 - minute breaks → ~7.5% CO₂

Despite heavy disturbance of the chamber atmosphere, **in-vessel conditions remained stable and close to the setpoint.**

Take-home message

- The CO₂ inside shaken vessels is remarkably **robust** against door openings.
- Even large chamber fluctuations cause only **small, slow changes** inside the tubes.
- Repeated openings produce a predictable equilibrium — In this study between 0.5–1.6% below the setpoint.
- This effect is highly **reproducible**, and with an appropriately higher CO₂ setpoint (via software control), the **deviation can be eliminated**.