

Kuhner shaker


FlowCon 2/3/4

Stand-alone gas mixing device



FlowCon 4

mixing up to 4 gases

Dimensions

FlowCon 2/3/4

The Kuhner FlowCon has been developed for the R&D environment and initial production steps in bioreactors with disposable bags. Whether stabilizing pH in cell cultivations with CO₂ or reducing oxygen concentration for microaerophilic organisms, the FlowCon gas mixer meets all user needs. The FlowCon can be used as a stand-alone device or can be integrated with the Kuhner equipment family (incubator shakers, disposable bioreactors).

Reproducibility

Up to four mass flow controllers ensure high accuracy and reproducibility of the gas mixture. The number of mass flow controllers can be adapted to individual requirements (FlowCon 2/3/4). An integrated pressure sensor prevents over pressure in the connected bioreactor by reducing the gas flow in the first step and shutting off the gas flow in the second step. Step limits can be freely programmed by the user. Integrated air filters and back pressure valves ensure interference free and long operating times for the FlowCon.

The Kuhner Insight software provides data logging, programming, controlling and calibration.

- + Mixing up to four gases (gas mixtures can also be connected)
- + Flow rates:
 - 0-1 [sL/min]
 - 0-2 [sL/min]
 - 0-20 [sL/min]
- + Integration in Kuhner shakers
- + Kuhner Insight Software control
- + Integrated pressure sensor

FlowCon 2/3/4

can be used for different types of bioreactors

Kuhner Insight Software

provides data logging, programming,
controlling & calibration



Technical data

<ul style="list-style-type: none"> Number of mass flow controllers 2, 3 or 4 	<ul style="list-style-type: none"> Stand-alone operation yes
<ul style="list-style-type: none"> Dimensions 430x225x158 (DxWxH) [mm] 	<ul style="list-style-type: none"> Software Kuhner Insight (optional): data logging, programming, controlling, calibration
<ul style="list-style-type: none"> Ambient conditions 5-40 [°C], max. 85 [% rH] 	<ul style="list-style-type: none"> Interface, standard CAN-Bus
<ul style="list-style-type: none"> Input up to 4 gases, push connector 6 [mm] + air + oxygen + nitrogen + carbon dioxide (also gas mixtures can be connected) 	<ul style="list-style-type: none"> Interface, optional USB, Ethernet, digital, analogue
<ul style="list-style-type: none"> Input pressure 1-2 bar 	<ul style="list-style-type: none"> Pressure sensor integrated, 0-400 [mbar], programmable, prevents over pressure in the connected vessel
<ul style="list-style-type: none"> Output 1, (tube fitting for tube ID 13 [mm]) 	<ul style="list-style-type: none"> Air filter integrated, prevents blocking of the mass flow controllers
<ul style="list-style-type: none"> Flow rates <ul style="list-style-type: none"> • 0-1 [sL/min] • 0-2 [sL/min] • 0-20 [sL/min] 	<ul style="list-style-type: none"> Back pressure valve integrated, prevents reflux into the mass flow controllers
<ul style="list-style-type: none"> Accuracy ±1.5% from the end value 	<p>Technical data subject to change</p>



Mixing up to 4 gases

www.kuhner.com

Kuhner

Represented by:

Adolf Kühner AG • since 1949

Dinkelbergstrasse 1
CH-4127 Birsfelden (Basel)
Switzerland

phone +41 (0) 61 319 93 93
fax +41 (0) 61 319 93 94
office@kuhner.com

www.kuhner.com