

Mammalian Cell Cultivation using AerationCaps for Shake Flask Gas Supply

Johannes Möller^{1,2,*}, Krathika Bhat¹, Tibor Anderlei^{3,+}, Christine Dambaur³, Ralf Pörtner^{1,2}



Objective

Mammalian cell cultivation using shake flasks requires agitation in CO₂ enriched atmosphere for pH control [1, 2]. To perform cultivations in rather simple shakers, e.g. only temperature and shaking control, AerationCaps (Kühner AG) are an add-on device to supply and maintain a controlled gas atmosphere needed for mammalian cell culture. Therefore, the applicability of AerationCaps is investigated in this poster for mammalian cell culture processes.



Aeration setup

A scheme of the Kühner AerationCap mounted on a shake flask is shown in Figure 1. The AerationCap has two luer-lock connectors to allow the gas supply. To prevent air from the surroundings entering into the flask, the AerationCap has a rubber sealing on the bottom. The caps can be easily installed on the shake flask and removed as needed, without any additional modifications to the flask. In conjunction with a controlled gas supply (i.e., FlowCon), the caps can be used to ensure the desired atmosphere inside the flasks, regardless of the incubator type.

To control the composition of the gas supply, Kühner AG offers FlowCon, a gas mixing device with up to four inputs and an adjustable output flow rate.

The output of the gas mixing-device can be equilibrated, humidified and divided between four flasks using the Kühner GasDivider. An exemplary set-up of the aeration system is shown in Figure 2.




Figure 2 Setup with AerationCap, GasDivider and FlowCon




Figure 1 Schematic of shake flask with AerationCap connected to gas supply by luer-lock attachments

1. FlowCon
2. GasDivider
3. AerationCap



Case Study: pH buffering in cell culture

CHO DP-12 cells were cultivated using AerationCaps (n=2) and FlowCon (flow rate: 120 ml min⁻¹) to maintain 5% CO₂ for pH buffering and humidity control. No aeration or humidity control was adjusted using the incubator controls. Cells from the same preculture were cultivated in a reference incubator with 5% CO₂ control (n=2). The Shake Flask Reader (SFR, PreSens) was used with the LT-X shaker (Kühner AG) for non-invasive, online monitoring of pH and dissolved oxygen (DO) in the shake flasks.

Comparison of cell growth using pH buffering through AerationCaps vs. pH buffering through incubator atmosphere




Figure 3 Cultivation of CHO DP-12 cells using aeration cap and FlowCon3 to maintain 5% CO₂ Flask volume: 80 ml, working volume: 40 ml, medium: TC-42 (Xell AG) with 6 mmol l⁻¹ glutamine




Figure 4 Experimental results of CHO DP-12 cultivations with and without AerationCaps; A: Cell growth, B: Glucose (solid line) and Lactate (dashed line) concentrations

✓ Using AerationCaps, **cell growth and metabolic activity are comparable** to the cell growth and metabolic activity at incubator atmosphere with 5% CO₂



Conclusion

FlowCon with AerationCaps can be used for controlled aeration of shake flasks and was successfully applied to pH buffering for cell culture

- ✓ Can be used in incubators without gas supply control as a **simple add-on**
- ✓ **Comparable pH and DO** maintained in the flasks with AerationCaps (seen in SFR data, not shown), **proven to support cell growth** similar to that in incubator-controlled gas supply

Contact

*johannes.moeller@tuhh.de, *tanderlei@kühner.com

¹Hamburg University of Technology, Institute of Bioprocess- and Biosystems Engineering Germany, 21073 Hamburg, Germany

²Cell Service Center, Hamburg University of Technology

³Adolf Kühner AG, 4217 Birsfelden, Switzerland

References

[1] Möller et al. (2019) DOI: 10.1007/s00449-019-02089-7

[2] Möller et al. (2020) DOI: 10.1016/j.compchemeng.2019.106693

Acknowledgement

Thanks to Prof. T. Noll (Bielefeld University) for providing the cell line used in this work. The authors would like to thank Kühner AG for providing the Shake Flask Reader, FlowCon and aeration caps, which enabled this study.