

LT-Z Peltier

Precise. Compact. Efficient.

NEW



Designed for your lab

Space is limited and valuable

Compact footprint. Stackable – also with other equipment. Double shaking capacity option.

Integration should be effortless

Seamless lab monitoring integration. Multi-bus compatible. Ready for diverse protocols.

Cleaning wastes time

Self-disinfecting Evaporation cup. Easy access. Easy handling. Autoclavable. Dishwasher-safe.

Silence matters

Whisper-quiet direct drive. High-grade wall insulation.



Efficient humidity control



Stackable unit

Perfect for your application

Reliability is non-negotiable

Shaking all day, every day. Maintenance-free direct drive.

Requirements change. The LT-Z adapts

Flexible shaking diameters. Smart accessory system. Upgradable with online technology.

Hassle-free validation. Ready for GMP use

Applied and tested IQ/OQ validation package. Customizable documents.

Hygiene is key

Clean water vapor humidification. Stainless steel interior. Design conforms to DIN standards.



GMP qualification



5 year warranty

New humidity control

Precise, reliable humidity control with Kuhner's patented **Cleva (Clean Evaporation) system**. Designed for GMP environments. Little maintenance. High level of hygiene.

- **Clean water vapor** – no mist, no condensation
- **Easy access** – tool-free removal
- **Easy handling** – magnetic positioning
- **GMP-compatible** – autoclavable & dishwasher-safe
- **Self-disinfecting** – no detergents required
- **Interchangeable** – clean one, use another



Upgrade your shaker with online monitoring


The new Kuhner microTOM elevates your screening on the fast track by measuring the growth kinetics of your cultures in each well of a 24 or 96 standard deepwell plate.

Kuhner microTOM



The Kuhner TOM, a fully automated setup, provides reliable and reproducible online data 24/7 without additional manual work.

Kuhner TOM



Want to learn more?
Contact us for more
information!

